

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06

Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://labconco.nt-rt.ru/> || bon@nt-rt.ru



Protector Pass-Through Laboratory Hoods

Protector Pass-Through Laboratory Hoods offer sashes on two sides allowing operation on either side of the fume hood. An audible alarm sounds to prevent simultaneous opening of both sashes. With unobstructed visibility from either side, this hood is ideal for conducting chemistry demonstrations. An optional right-side or left-side viewing window further enhances visibility. These hoods may also be installed in the wall between two rooms for sample preparation.

Use Our Product Selection Tools

Scout™ Lab Equipment Selector

Let Scout lead you to the right laboratory hood for your application.

Features

- By-pass airflow design and by-pass blocks, one for each sash
- Glacier white powder-coated steel and aluminum exterior
- Powder-coated steel Eco-Foil™ air foil with Clean-Sweep™ airflow openings.
- Cord-Keeper™ slots on left and right side of air foil
- Chemical-resistant molded composite panel interior sides and pre-set baffle with flame spread less than 25 per ASTM E84
- Powder-coated aluminum sash handle with Clean-Sweep™ openings and Clean-Sweep™ slots on the powder-coated steel sash tracks
- Two tempered safety glass vertical-rising sashes with cable pulleys and sash stops mounted at 25" (63.5 cm) opening height
- 37.5" (95.3 cm) high sightline from the work surface to the header panel
- Audible alarm that sounds when both sashes are open at the same time
- Removable front and side panels and front access panels for access to plumbing and electrical wiring
- Pre-wired T8 fluorescent lighting with vapor-proof design, one on each side and ADA-compliant light and blower switches on the front side
- Factory prepared for up to eight service fixtures, four on each side
- Powder-coated stainless steel 12.81" ID exhaust connection

Side Window Models Feature

- Right or left side-mounted tempered safety glass window, 14.9" x 23.9" (37.8 x 60.7 cm)

Fixtured Models May Features

- Two pre-plumbed service fixtures with forged brass valves with brass tubing for gas and copper tubing for cold water. Components for converting either or both fixtures to air or vacuum are provided. Inlet tubing is not provided.
- One pre-wired GFCI electrical duplex receptacle on lower right side

Required Accessories (Not Included)

- Remote blower
- Ductwork
- Work surface
- Base cabinet or stand

Optional Accessories

- Service Fixture Kits (some models include service fixtures)
- Electrical Duplex Kits (some models include electrical duplex)
- Guardian Airflow Monitor Kits
- Sash Stop Kits

4' Protector Pass-Through Laboratory Hood with 2 fixtures, 1 receptacle, 115V



Catalog # 113400002

4' Protector Pass-Through Laboratory Hood with 2 fixtures, 1 receptacle, left window, 115V



Catalog # 113407002

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Electrical Duplexes: 1

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
1 receptacle, right window, 115V**



Catalog # 113408002

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Electrical Duplexes: 1

Service Fixtures: 2

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
230V**



Catalog # 113400021

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Service Fixtures: 2

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
left window, 230V**



Catalog # 113407021

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Electrical Duplexes: 1

Service Fixtures: 2

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
115V**



Catalog # 113400001

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
left window, 115V**



Catalog # 113407001

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
right window, 115V**

Catalog # 113408001

Nominal Width: 4'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Service Fixtures: 2
Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood with 2 fixtures,
right window, 230V**



Catalog # 113408021

Nominal Width: 4'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Service Fixtures: 2
Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood, 230V**



Catalog # 113400020

Nominal Width: 4'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood, left window,
230V**



Catalog # 113407020



Nominal Width: 4'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood, 115V**



Catalog # 113400000

Nominal Width: 4'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood, left window,
115V**



Catalog # 113407000

Nominal Width: 4'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood, right window,
115V**



Catalog # 113408000

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Style: Benchtop

**4' Protector Pass-Through
Laboratory Hood, right window,
230V**



Catalog # 113408020

Nominal Width: 4'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
1 receptacle, 115V**



Catalog # 113500002

Nominal Width: 5'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Electrical Duplexes: 1

Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
1 receptacle, left window, 115V**



Catalog # 113507002

Nominal Width: 5'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Electrical Duplexes: 1

Service Fixtures: 2

Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
1 receptacle, right window, 115V**



Catalog # 113508002

Nominal Width: 5'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required

Electrical Duplexes: 1

Service Fixtures: 2

Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
115V**



Catalog # 113500001

Nominal Width: 5'

Region: U.S. and Canada,
International

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
230V**



Catalog # 113500021

Nominal Width: 5'

Region: U.S. and Canada,
International

Blower Requirements: Remote
blower required
Style: Benchtop

International
Blower Requirements: Remote
blower required
Service Fixtures: 2
Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
left window, 115V**



Catalog # 113507001

Nominal Width: 5'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
left window, 230V**



Catalog # 113507021

Nominal Width: 5'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Service Fixtures: 2
Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
right window, 115V**



Catalog # 113508001

Nominal Width: 5'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Style: Benchtop

**5' Protector Pass-Through
Laboratory Hood with 2 fixtures,
right window, 230V**



Catalog # 113508021

Nominal Width: 5'
Region: U.S. and Canada,
International
Blower Requirements: Remote
blower required
Service Fixtures: 2
Style: Benchtop

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06

Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://labconco.nt-rt.ru/> || bon@nt-rt.ru