

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06

Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://labconco.nt-rt.ru/> || bon@nt-rt.ru

CAapture BT Benchtop CA Fuming Systems



The CAapture™ BT is a spacious Cyanoacrylate (CA) evidence fuming chamber with a durable easy-to-clean stainless steel liner. It is designed with input from Forensics Experts to protect you and your evidence.

- Programmable operating system automates process, allowing control over critical parameters
- Acti-Vent Smart Controls protect evidence from over-fuming and protect personnel
- Eye-level LCD display allows precise parameter reading
- Pre-Humidity Incubation and fuming temperature up to 424°F
- Dual blower design ensures airflow patterns that optimize fuming and exhaust cycles

Spacious Cyanacrylate (CA) evidence fuming chambers for the forensic lab benchtop or portable stand.

Additional Features & Specifications

This benchtop chamber with microprocessor control automates the Cyanoacrylate fuming of evidence. The

monitoring of all critical parameters in the fuming process with the precision to store and easily repeat specified programs. Its unique interior space with shelving system maximizes fuming capacity and accommodates large evidence including long guns.

Features & Benefits

- 27" w x 27" d x 44.5" h interior.
- Durable, easy-to-clean type 304 brushed stainless steel interior, epoxy-coated steel exterior and tempered safety glass door.
- Electronic door lock with power outage manual override prevents opening during the fuming cycle.
- Interior shelving support accommodates up to four shelves, eight half shelves or 16 rods.
- Dependable, energy-efficient ECM recirculation blower moves air through the chamber during the humidifying and fuming cycles.
- Exhaust blower removes air from the chamber through slots located on the back, right side of the cabinet after passing through a pre-filter and carbon filter during the exhaust or purge cycle.
- Full-color LCD control screen is located at eye level for at-a-glance process monitoring. Displays number of completed fuming cycles, filter life remaining, low-level water reservoir alerts and open door alerts.
- Microprocessor-controlled operating system provides flexibility and control:
- Adjustable chamber humidity to 80% +/- 5% relative humidity
- Adjustable pre-fuming humidity incubation for old evidence hydration
- Adjustable fuming temperature up to 424° F within +/-1%
- Adjustable fuming and exhaust times
- Stores up to 20 programs

- Standby exhaust mode to refresh chamber air when not in use
- Diagnostic mode to calibrate sensors
- Audible/visual alert for program cycle completion
- Acti-Vent™ Smart Controls purge air to protect evidence and the user. Purging continues up to 30 minutes after exhaust cycle ends. Purging may be programmed to continue during evidence unloading. Cleaning mode enables air to be exhausted when the door is open to protect the user from exposure to noxious vapors.
- Convenient drawer below the chamber, 28" w x 13" d x 1.5" h, has flushmounted disappearing handle, platform for supporting tins while adding CA, and compartments to store glue, hanging rods and other supplies.
- Includes four stainless steel hanging rods, interior LED light, exhaust pre-filter, exhaust carbon filter, four leveling feet, power cord and plug.
- Optional accessory CApture BT Mobile Stand adds portability and storage.

[CApture BT Chamber, 115V North America](#)



Catalog # 3170000

Product Type: CA Fuming System
Nominal Width: 3'
Region: U.S. and Canada
Power Cord & Plug: North America, 115 volts

[CApture BT Chamber, 230V North America](#)



Catalog # 3170001

Product Type: CA Fuming System
Nominal Width: 3'
Region: U.S. and Canada
Power Cord & Plug: North America, 230 volts

[CApture BT Chamber, 230V Schuko](#)



Catalog # 3170002

Product Type: CA Fuming System
Nominal Width: 3'
Region: International
Power Cord & Plug: Schuko

[CApture BT Chamber, 230V British \(UK\)](#)



Catalog # 3170003

Product Type: CA Fuming System
Nominal Width: 3'
Region: International
Power Cord & Plug: British (UK)

CApture BT Chamber, 230V
China/Australia



Catalog # 3170004

Product Type: CA Fuming System
Nominal Width: 3'
Region: International
Power Cord &
Plug: China/Australia

CApture BT Chamber, 230V India



Catalog # 3170005

Product Type: CA Fuming System
Nominal Width: 3'
Region: International
Power Cord & Plug: India

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06

Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://labconco.nt-rt.ru/> || bon@nt-rt.ru