

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06

Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

Киргизия (996)312-96-26-47

Россия (495)268-04-70

Казахстан (772)734-952-31

<https://labconco.nt-rt.ru/> || bon@nt-rt.ru

ScrubAir Pipette Washer/Dryers

The ScrubAir Pipette Washer/Dryer [saves 98% on water and takes 85% less time to use.](#)



Compared to traditional siphoning washing, rinsing and oven drying, the ScrubAir consumes less water, time and money, requiring as little as 3 1/4 hours to wash, rinse and dry up to 60 pipettes, compared to 23 hours or more with the traditional method. Eco friendly, the ScrubAir Standard Wash Cycle uses only 12.5 L of tap or pure water, compared to 600 L used with a siphoning method. Three cycle programs are factory-set but may be user-modified. Start time may be delayed up to 8 hours. The LCD display shows both the cycle segment and total time remaining.

- Saves water, time and money
- Automated, requires far less glassware handling than traditional methods
- Fast and thorough cleaning with percolating action
- Programmable cycle settings including delayed start

The ScrubAir Pipette Washer/Dryer saves 98% on

pipette washing methods.

Additional Features & Specifications

The ScrubAir Pipette Washer/Dryer uses water, detergent and compressed air to wash, rinse and dry volumetric pipettes. Water and detergent are directly injected into the pipettes, then compressed air is introduced into the water line producing a percolating action that scrubs the interior walls. The water is then propelled out the top of the pipettes, washing the exterior of the pipettes as it cascades down. After washing and rinsing, the ScrubAir dries the pipettes in two to eight hours. No separate oven is required.

This automatic washer/dryer can accommodate pipettes ranging in size from 0.135" (3.4 mm) stem diameter and 7.5" long up to 0.34" (8.6 mm) stem diameter and 25" (635 mm) long. Depending on the bulb and stem diameters, up to 60 small or 12 large pipettes may be washed at one time. Heated models also include a built-in 750 watt heating element, activated during the wash, rinse and dry cycles, and a pure water inlet to select up to 10 pure water rinses.

All models include brushed Type 304 stainless steel and powder-coated steel exterior, Type 304 stainless steel pipette holder, Type 304 stainless steel back splash, clear acrylic cover, feedwater inlet/drain hose with quick connect hose coupling and 2 aerator adapters for connection to standard laboratory faucets, rear-mounted polyethylene tube for attachment to regulated air source or other suitable inert gas, and 1/8" ounce dispense pump included for attachment to liquid detergent container (Detergent sold separately). Models with built-in heater also include rear-mounted polyethylene tube for attachment to pure water source

All models require tap/pure feedwater with 3/8" OD (minimum) water line, regulated compressed air with 3/8" OD air line, drain plumbing 1.25" (3.2 cm) ID or larger, and liquid detergent (LabSolutions Liquid Detergent recommended). See brochure for water temperature, pressure, and airflow pressure requirements.

[ScrubAir Pipette Washer/Dryer with Heater 115V US](#)



Catalog # 4612021

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada, International

[ScrubAir Pipette Washer/Dryer with Heater 230V US](#)



Catalog # 4612041

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada, International

[ScrubAir Pipette Washer/Dryer with Heater 230V Schuko](#)



Catalog # 4612031

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada, International

[ScrubAir Pipette Washer/Dryer with Heater 230V UK](#)



Catalog # 4612061

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada, International

[ScrubAir Pipette Washer/Dryer with Heater 230V AUS/CHN](#)



Catalog # 4612071

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada, International

[ScrubAir Pipette Washer/Dryer with Heater 230V India](#)



Catalog # 4612081

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada, International

[ScrubAir Pipette Washer/Dryer](#)
[115V US](#)



Catalog # 4612020

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada,
International

[ScrubAir Pipette Washer/Dryer](#)
[230V US](#)



Catalog # 4612040

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada,
International

[ScrubAir Pipette Washer/Dryer](#)
[230V Schuko](#)



Catalog # 4612030

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada,
International

[ScrubAir Pipette Washer/Dryer](#)
[230V UK](#)



Catalog # 4612060

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada,
International

[ScrubAir Pipette Washer/Dryer](#)
[230V AUS/CHN](#)



Catalog # 4612070

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada,
International

[ScrubAir Pipette Washer/Dryer](#)
[230V India](#)



Catalog # 4612080

Product Subcategory: ScrubAir
Style: Benchtop
Region: U.S. and Canada,
International

Архангельск (8182)63-90-72
Астана (7172)727-132
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89
Иваново (4932)77-34-06

Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Липецк (4742)52-20-81

Киргизия (996)312-96-26-47

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16

Россия (495)268-04-70

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13

Казахстан (772)734-952-31

Сургут (3462)77-98-35
Тверь (4822)63-31-35
Томск (3822)98-41-53
Тула (4872)74-02-29
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Ярославль (4852)69-52-93

<https://labconco.nt-rt.ru/> || bon@nt-rt.ru